

*iLase\**



The Powerful, Affordable, Pocket Portable  
Addition to Your Practice

**BIOLASE**

WITH NO FOOT PEDAL, POWER CORD  
OR EXTERNAL CONTROLS,



THE REVOLUTIONARY iLASE™ TAKES DENTAL LASERS TO  
**POCKET PORTABLE**



"I use the iLase daily. Because of the small size, it is always ready to go. I use it for cleaning extraction sockets, uncovering implants and sulcus contouring around implants (prevents cement from getting on the implant), plus much more!"

*– Dr. Bret Dyer, Sugarland, Texas*

## The First Personal Laser

All the power and control is concentrated in your hand for uniquely personal performance.

- ✦ No foot pedals to clutter your operatory floor
- ✦ Avoid getting tangled in cords and wires with single, self-contained handpiece
- ✦ Increase productivity with an unrestricted, truly personal experience

## Versatility at the touch of a finger

With the ability to engage the finger switch from different positions, the iLase allows maximum comfort and working efficiency, whether performing procedures in the anterior or posterior of the mouth.

- ✦ Wrap around finger switch accessible from any angle in the oral cavity
- ✦ "Standby" & "Ready Mode" protection
- ✦ Fully autoclavable for easy infection control

## Clinically Convenient

The iLase amazes because so much convenience is packed into such a small, portable device. Surgical tips are disposable for easy setup and teardown. And the iLase ships with 10 factory-loaded surgical presets, with slots for 14 more of your favorite procedures.

### PRESETS INCLUDE:

- ✦ Gingival troughing for crown impressions
- ✦ Gingivectomies
- ✦ Treatment of perio pockets
- ✦ Sulcular debridement
- ✦ Frenectomies
- ✦ Excisions
- ✦ Aphthous ulcers
- ✦ Hemostasis
- ✦ Exposure of unerupted teeth
- ✦ Implant recovery
- ✦ Crown lengthening

### TIP FEATURES:

- ✦ Disposable single-use tips
- ✦ Bendable & adjustable
- ✦ Full assortment of lengths & diameters for all clinical situations
- ✦ Pre-initiated tips available



## Balanced and Ergonomic

The precise ergonomics of the iLase handpiece also reduce arm, wrist, and hand fatigue for clinicians.

## Power That Lasts

The iLase delivers power when you need it most – at the time of treatment. Whether the procedure is long or short, iLase has the battery life to meet your needs.

- ✦ Delivers up to 5 Watts of peak power for even the toughest soft-tissue cases
- ✦ Adjustable power in tenths of a Watt for clinical precision
- ✦ Convenient countertop battery charger and extra batteries included
- ✦ 3 Watts of continuous power



The iLase charging station recharges up to four iLase batteries simultaneously.

## Made in the USA

Unlike most medical and dental lasers, BIOLASE proudly designs and manufactures all of our lasers in the USA at our facility in Irvine, California.

## Affordable for All

The iLase was designed to be affordable, so that each dentist and hygienist can work with his or her own laser. One more way the iLase is taking truly personal performance to a whole new level.



## Training & Certification

Online training by experienced diode laser dentists is available through BIOLASE in collaboration with the WCLJ, the world's largest dental technology education organization. To learn more about the WCLJ, or to schedule training, visit [www.wcli.org](http://www.wcli.org)

# *ilase\**



## BIOLASE

For more information, visit [biolase.com](http://biolase.com)   

### BIOLASE HEADQUARTERS

4 Cromwell  
Irvine, CA 92618, USA  
Toll-Free: (888) 424-6527  
Phone: (949) 361-1200  
Service: (800) 321-6717  
Fax: (949) 273-6687

### EUROPE

BIOLASE Europe GmbH  
Paintweg 10  
92685 Floss, Germany  
Phone: +499603808252  
Fax: +499603808250

### INDIA

BIOLASE Technology India Pvt. Ltd.  
404, Agarwal Golden Chambers  
Fun Republic Road, Off New Link Road  
Andheri (W), Mumbai — 400 053  
Phone: +91 22 6195 3400